



Case study: Empowering Iga to enjoy daily life

Young Iga's lack of muscular tone means she struggles to maintain a position without losing trunk control.

Now she can sit longer, and more independently. As a result, she can participate in and focus on daily activities, with more autonomy. Her fine motor skills are improving too: she's now able to eat on her own!

Iga's improvement is down to a combination of our innovative vacuum posture cushions.

Central to her improved seated position and trunk control is Base Doro.

The cushion's prominent pommel front wedge places the legs in a neutral position preventing contracture and encouraging work to maintain a proper torso position. The high lateral sides provide additional trunk support.

For Iga, we addressed a posterior pelvic tilt by correct shaping of the granules to guide and stabilise her pelvis into a more neutral position. We then used our slim Armrest cushion to create additional lumbar support and our Wedge to support her feet, even when Iga's sat in a standard chair.

You can find out more about Base Doro here: www.aatgb.com/base-doro-2/

There is more information about our vacuum posture cushions here: www.aatgb.com/posture-cushions/

Iga before...



Iga after...



AAT GB Ltd

www.aatgb.com

e: sales@aatgb.com

t 01327 340083